

# Breastmilk Feeding & Storage Solutions




## Breastmilk Bottle Sets

Safe, BPA-free plastic or glass breastmilk bottles are designed for long-term breastmilk storage and retain breastmilk's beneficial properties.

- Compatible with Medela breastpumps\* to pump, store and feed with one bottle.
- New wide base silicone nipples.
- Available in convenient sizes - 5 oz & 8 oz polypropylene and 4 oz & 8 oz glass.
- Volume marks in ounces and milliliters.
- Dishwasher and microwave safe.
- Glass breastmilk bottles are recyclable.



## Wide Base Nipples



- Natural feel and shape.
- Available in slow and medium flow for feeding newborns and infants.
- Compatible with all Medela breastmilk bottles\*\*.
- Dishwasher and microwave safe.



## Pump & Save™ Breastmilk Bags

Available in 20-pack & 50-pack.

- Unique design – pump directly into bag for fast, easy breastmilk storage.
- Compatible with Medela breastpumps\*.
- Double-walled for safe and long-term breastmilk storage.
- No-leak, easy-to-close zipper top.
- No-spill, stand-up bottom.



## Breastmilk Freezing & Storage Bulk Pack

(12) 2.7 oz/80 mL bottles per pack.

- Perfect for smaller amounts of pumped breastmilk.
- Easy storage in refrigerator or freezer.
- Compatible with Medela breastpumps\* for pumping directly into the container.



## Breastmilk Feeding Gift Set™

Value Pack

- Everything you need to collect, store and feed breastmilk.
- Compatible with Medela breastpumps\*.
- Includes: (3) 5 oz breastmilk bottles with slow flow nipples, collars & caps, (6) 2.7 oz/80 mL bottles and (20) Pump & Save breastmilk bags.



## Breastmilk Bottle Spare Parts

Always have clean parts on hand.

- Convenient mesh bag for cleaning bottle parts in the dishwasher – and keeping them together!
- For use with Medela wide base nipples and breastmilk bottles.
- Includes three sets of travel caps, collars, discs and solid lids - perfect for long-term breastmilk storage.\*\*

\*\*Requires wide base collar. Call Medela at 1-800-435-8316 if you need a wide base collar.

**Make sure breastfeeding is well established before introducing a nipple.**

\*2.7 oz/80 mL breastmilk bottles and 8 oz glass breastmilk bottles not for use with the Single Deluxe™ breastpump. Pump & Save breastmilk bags and 8 oz glass breastmilk bottles not for use with the Freestyle® breastpump.

Medela, Inc., P.O. Box 660, 1101 Corporate Drive, McHenry, IL 60051-0660 Phone: (800) 435-8316 or (815) 363-1166  
Fax: (815) 363-1246 Email: [customer.service@medela.com](mailto:customer.service@medela.com) [www.medela.com](http://www.medela.com)

Pump & Save, Pump-Store-Feed, Single Deluxe and Breastmilk Feeding Gift Set are trademarks and Medela and Freestyle are registered trademarks of Medela, Inc.

1908101 B 0310 © 2010 Medela, Inc. Printed in the USA.





# Breastmilk Feeding & Storage Solutions

Everything you need for safe breastmilk collection, storage and feeding.



## Breastmilk Feeding & Storage Solutions

If you are planning to store pumped breastmilk to feed your baby, Medela offers a complete line of BPA-free products that are designed for safe feeding and long-term storage of breastmilk.

Medela breastmilk bottles and Pump & Save™ bags are compatible with Medela breastpumps\* to pump, store and feed breastmilk with the same container. Use the chart below to help you select the perfect feeding and storage products for your situation.



medela	Breastmilk Freezing & Storage Containers			Breastmilk Bottles		Pump & Save Breastmilk Bags*	
	2.7 oz/80 mL			5 oz & 8 oz Polypropylene or 4 oz & 8 oz Glass*			
Long Term Storage (freezer)	■			■		■	
Pump • Store • Feed™ with one container	■ <small>(nipples sold separately)</small>			■			
Pump • Store						■	
Complete breastmilk feeding with nipple, travel cap and lid				■			
Feed baby breastmilk	■			■			
Durable and recyclable	■			■			
Disposable						■	
BPA-Free	■			■		■	

Freshly Expressed Breastmilk Storage Guidelines <i>(For Healthy Term Babies)</i>				
Room Temperature	Cooler with 3 Frozen Ice Packs	Refrigerator	Freezer	Thawed Breastmilk
4–6 hours at 66–78 °F (19–26 °C)	24 hours at 59 °F (15 °C)	3–8 days at 39 °F or lower (4 °C)	6–12 months 0–4 °F (-18–-20 °C)	use within 24 hrs

References: [www.BreastmilkGuidelines.com](http://www.BreastmilkGuidelines.com)

If your baby was born prematurely, these guidelines may differ slightly. You should check with your health care provider for the recommended storage guidelines for your specific situation. Never microwave breastmilk. Microwaving can cause severe burns to baby's mouth from hot spots that develop in the milk during microwaving. Microwaving can also change the composition of breastmilk.

For more information, or to find a lactation consultant near you, call our Breastfeeding National Network (BNN) at 1-800-TELL YOU or visit [www.medela.com](http://www.medela.com). Make sure breastfeeding is well established before introducing a nipple.

\*2.7 oz/80 mL breastmilk bottles and 8 oz glass breastmilk bottles not for use with the Single Deluxe™ breastpump. Pump & Save breastmilk bags and 8 oz glass breastmilk bottles not for use with the Freestyle® breastpump.



Medela, Inc.  
P.O. Box 660  
1101 Corporate Drive  
McHenry, IL 60051-0660